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ENGLISH / US

Sterile: Contents sterile unless package is opened, damaged or expired.

- A. INTENDED PURPOSE:** The Ranfac Bone Marrow Biopsy Needle is intended for the purpose of harvesting bone marrow specimens.
- B. INTENDED USER:** The Ranfac Bone Marrow Biopsy Needle is intended to be used by a licensed physician familiar with the possible side effects, typical findings, limitations, indications, and contraindications of bone marrow biopsy.
- C. INTENDED PATIENT POPULATION:** The Ranfac Bone Marrow Biopsy Needle is intended to be used on any patients that are suitable for bone marrow harvesting as determined by the physician.
- D. INDICATION(S):** The Ranfac Bone Marrow Biopsy Needle is indicated for the surgical harvesting of bone marrow.
- E. CONTRAINDICATION(S):** The Ranfac Bone Marrow Biopsy Needle should not be used for sternal procedures.
- F. INTENDED LIFETIME:** The Ranfac Bone Marrow Biopsy Needle is a sterile, single-use device intended for
 - Transient use (<60 minutes)
- G. INTENDED CLINICAL BENEFIT:** The Ranfac Bone Marrow Biopsy Needle is intended to harvest bone marrow specimens in order to minimize the need for external manipulations of the needle, to help recover adequate quantity and better quality of specimen (lower risk of sample loss) and to minimize the likelihood of complications such as patient discomfort, pain, or injury.

WARNINGS/CAUTIONS: For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. Do not use if package is open or damaged and if the expiry date has been exceeded. Do not continue to use if the device gets damaged during procedure. Do not bend the needle excessively, it may result to damage/breakage of the device. After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

Note: These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

H. DIRECTIONS FOR USE

1. Inspect the package for proper integrity. If undamaged, open the package using aseptic technique.
2. Remove the contents from the packaging.
3. Remove needle guard from needle and inspect the cutting edges of both cannula and stylet for any damage or other imperfections that would prevent the proper operation of the needle.
4. Following proper placing technique, locate the needle adjacent to the biopsy site.
5. Using gentle, but firm, pressure advance the needle by rotating it in an alternating clockwise - counterclockwise motion. Entrance into the marrow cavity is generally detected by decreased resistance.
6. Unlock the stylet cap by rotating it counterclockwise and remove the stylet from the cannula.
7. If included use the needle vise provided to cap the stylet after use for safe disposal.

If blood begins to flow out through the proximal end of needle and a luer cap is included, follow steps 8-9. If not, skip to step 10 or 12.

8. Install luer cap onto luer threads at proximal end of needle.
9. Once you are ready to continue with the procedure, remove cap.

Note: If cap is on too tight to be removed by hand, a flat head screw driver can be used.

If both aspirate and a trephine biopsy are required, it is recommended that separate needles be used in order to preserve sample quality. However, if using the needle to acquire bone marrow aspirate, follow steps 10-11.

10. Attach a male luer-lock syringe to the needle using the luer lock connection.
Note: Avoid using glass syringes because they may result in an inadequate seal and may break or fragment.
11. Apply negative pressure by quickly withdrawing the syringe plunger. Disengage the syringe and remove the aspirated sample.

CAUTION: Avoid excessively deflecting needle during use as this may result in the inability to acquire bone marrow aspirate.

If using needle to acquire a bone marrow sample, continue from step 12.

12. Very slowly and gently, advance the needle using clockwise-counter-clockwise motion until adequate marrow is obtained.

13. Insert the marked probe without force to check sample length in needle lumen. Remove the probe.

If needle handle includes a grey lever, follow steps 14-17 to acquire bone marrow sample.

14. Deploy the gray lever by moving it from the "OPEN" to the "CLOSED" position.
15. Slowly withdraw the needle from the patient.
16. Rotate the gray lever from the "CLOSED" to the "OPEN" position.
17. Place the probe guide onto the distal tip of the cannula. Insert the probe into the distal end of the probe guide. Advance the probe into the needle and collect the biopsy specimen as it emerges from the opening in the handle (proximal end).

CAUTION: If any difficulty is encountered during the transfer of the sample, DO NOT force the sample. Remove the probe, activate the grey lever to the "CLOSED" position and return it to the "OPEN" position. Re-insert the probe and transfer the specimen out the proximal end of the needle.

If an accessory is provided, follow steps 18-20 to acquire bone marrow sample.

18. Insert accessory fully into needle cannula.
19. Slowly rotate and remove the accessory and needle cannula together with the same clockwise-counterclockwise motion.
20. Remove the accessory from the needle cannula and push the specimen out with the probe.

If the product includes neither a grey lever nor an accessory, follow steps 11-24 to acquire bone marrow sample.

21. With a rapid motion completely rotate the needle ten times to the right and ten times to the left while applying light downward pressure to allow the tapered tip to help sever the biopsy specimen.
22. Pull the needle back two or three millimeters and with minimal pressure, direct its tip at a slightly different angle. Advance the needle two or three millimeters further, and completely rotate to the right two more times. This technique further assists to sever the specimen before withdrawing the needle.

CAUTION: Applying too much pressure when redirecting the needle may bend the needle.

23. Slowly remove the needle with the same clockwise-counterclockwise motion.
24. Place the probe guide onto the distal tip of the cannula. Insert the probe into the distal end of the probe guide. Advance the probe into the needle and collect the biopsy specimen as it emerges from the opening in the handle.

If needle vise is included, use step 25.

25. Use the needle vise provided to cap the needle after use for safe disposal.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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MedEnvoy Global B.V.